



Migraine & Light Sensitivity Eyewear

Semifinished Pucks

Can be surfaced, edged, and mounted into custom frames.

Materials: Polycarbonate, Trivex 1.6

Base Curves: 2, 4, 6, 8

- Avulux Rx Pucks Single Vision or Progressive

Materials: Polycarbonate, Trivex 1.6

Base Curve: 2, 4, 6, 8

- Avulux Plano Pucks

Materials: Polycarbonate

Base Curve: 6

Finished Glasses

As seen on Avulux.com

- Avulux Acetate Frames with Plano Avulux Lenses
- Avulux Acetate Frames + Avulux Titanium Clip-on containing Plano Avulux Lenses (clear Rx lenses may be surfaced, edged, and mounted into the Acetate frame)

Where to Order

Avulux Rx Semifinished/Plano Lenses | contact sales@avulux.com

Avulux Rx Surfaced Lenses | refer to Authorized Avulux Lab list
<https://avulux.com/pages/ecp-authorized-lab>

Finished Avulux Branded Glasses | contact sales@avulux.com